Flaem Original Accessories for Nebulizers



We recommend the use of FLAEM genuine accessories, carefully designed and manufactured so to assure therapy effectiveness and safety.

Soft-Touch Flaem Masks

Safe & comfortable with perfect face adherence, the Soft-Touch masks, have the external edge made in "TPE", a very soft and fully biocompatible material. The anti-dispersion system is optimizing the therapy, by reducing the drug dispersion during the inspiratory and expiratory phases.

						Neb					
	Code	Model	Colour	Package type	RF7 Plus	RF7 Dual Speed	RF6 Plus	RF6 Basic	RF8 Plus	RF8 Basic	Special pack available
-0	ACO611P	Mask soft Adult + elastic band	•	Single pack	*	*	~	~	*	~	50 pieces package (code: ACO473P)
🕲	ACO612P	Mask soft Child + elastic band	•	Single pack	*	*	~	~	~	~	50 pieces package (code: ACO472P)
-	ACO437P	Mask soft Adult + elastic band	•	Single pack	~	~	~	~	~	~	50 pieces package (code: ACO483P)
	ACO438P	Mask soft Child + elastic band	•	Single pack	~	*	~	~	~	~	50 pieces package (code: ACO482P)
	ACO462P	Mask soft Adult + Child + elastic band	•	Single pack	~	~	~	~	~	~	-

Masks Soft-Touch made in Flaem - Baby

Full range of soft & safe silicone masks for babies. The masks can be connected with the universal L-bend connector, that is allowing to reach the baby in every position.

			-		Nebulizer compatibility						
	Code	Model	Colour	Package type	RF7 Plus	RF7 Dual Speed	RF6 Plus	RF6 Basic	RF8 Plus	RF8 Basic	
01	ACO389P	FLAEM BABY Mask Kit (Mask size 0 with bend) Size 0 = Neonatal - 1 month	Silicon	Single pack with manual	~	~	~	~	~	~	
01	ACO390P	FLAEM BABY Mask Kit (Mask size 1 with bend) Size 1 = 1-11 months	Silicon	Single pack with manual	~	*	~	~	~	~	
	ACO391P	FLAEM BABY Mask Kit (Mask size 2 with bend) Size 2 = 1-3 years	Silicon	Single pack with manual	~	*	~	~	~	~	
	ACO013B	L-Bend BABY connector	-	50 pieces	~	~	~	~	~	~	

Mouthpiece	es										
Flaem original mouthpieces	Code	Model	Colour	Package type	RF7 Plus	RF7 Dual Speed	She Plus	ompatik KŁę Basic R	bility Snld 848	RF8 Basic	
	ACO441P	Mouthpiece with expiratory valves	•	50 pieces	~	-	~	-	~	-	
	ACO427P	Mouthpiece	-	50 pieces	-	~	-	~	-	~	
	ACO250P	Mouthpiece+Nosepiece	•	1	~	-	~	-	~	-	



Non-Invasive with valves



Nosepieces

Flaem original mouthpieces	Code	Model	Colour	Package type	RF7 Plus	RF7 Dual Speed	RF6 Plus	RF6 Basic	RF8 Plus	RF8 Basic	
e e	ACO448P	Nosepiece Non-Invasive ADULT/PEDIATRIC	•	50	~	-	~	-	~	-	
	ACO450P	Nosepiece Standard ADULT	-	50	~	~	~	~	~	~	
Ni li	ACO451P	Nosepiece Standard PEDIATRIC	-	50	~	~	~	~	~	~	

Filters - Tubing - Controls

					Neb					
Flaem original filters, tubing and controls	Code	Model	Package type	RF7 Plus	RF7 Dual Speed	RF6 Plus	RF6 Basic	RF8 Plus	RF8 Basic	Special pack available
	ACO428P	Drug Saver Manual control	50	~	~	~	~	~	~	-
\bigcirc	ACO225P	Nebulizer tube 2 mt. with fittings	1	~	~	~	~	~	~	25 pieces (code: ACO457P)
	ACO35	Nebulizer tube 1 mt. with fittings	1	~	~	~	~	~	~	50 pieces (code: ACO432P)
\bigcirc	ACO578P	Nebulizer tube 2mt. with Long-therapy fittings. With TPU fitting.	1	~	~	~	~	~	~	-
l.	ACO164P	Kit 2 Filters POREX Professional (type XS-9226	2	-	-	-	-	-	-	50 pieces (code ACO455P)
150	ACO439P	Kit 2 Filters POREX Basic PE- PET	2	-	-	-	-	-	-	50 pieces (code: ACO485P)
<u>J</u> e	ACO45	Kit 2 Filters POREX Professional (type X-6101)	2	-	-	-	_	-	-	50 pieces (code: ACO456P)