

	Code	Model	Details / Kit composition	Suction Unit compatibility				
				Aspira-Go 30lt	Aspira-Go 30lt	Aspira 20lt	Apsira 13.5	Aspira 11lt
	ACO002A	Autoclavable JAR 1000 ml	Reusable JAR in polycarbonate silicon-sealed cover and antireflux valve	✓	✓	✓	✓	✓
	ACO004A	Set with autoclavable JAR, 1000ml filter and tubing	Reusable JAR in polycarbonate 1000ml, filters, tubes	✓	✓	✓	✓	✓
	ACO001A	Single patient JAR 800ml	Single-patient JAR with silicon-sealed cover and overflow valve	✓	✓	✓	✓	✓
	15224A0	Additional Jar-Holder	Universal left or right Jar Holder, complete with antivibrating feet	✓	✓	✓	✓	✓
	ACO555P	Set of 2 filters	Antibacterial, antiviral, hydrophobic filters (single-patient)	✓	✓	✓	✓	✓
	ACO554P	Set 3 silicon tubes	Autoclavable set of 3 tube	✓	✓	✓	✓	✓
	ACO003A	Set 2 silicon tubes	Autoclavable set of 2 tube	✓	✓	✓	✓	✓
	18772	Manual vacuum flow control	Disposable manual suction regulator	✓	✓	✓	✓	✓
	16536B0	Protective carrying case	Carrying case for ASPIRA-GO models, with reflecting inserts	✓	✓	-	-	-
	ACO653P	Li-ion battery set	Replaceable Li-ion battery set for ASPIRA-GO models, fast to charge and long runtime	✓	✓	-	-	-
	ACO203P	Universal switching power adapter	Universal switching power adapter with EU plug (100-240V 50/60Hz)	✓	✓	-	-	-