






# Accessories for Nebulizers



## Soft-Touch Flaem Masks

Safe & comfortable with perfect face adherence, the Soft-Touch Masks, have the external edge made in "TPE", a very soft and fully biocompatible material. The anti-dispersion system is optimizing the therapy, by reducing the drug dispersion during the inspiratory and expiratory phases.




Code	Description	Package type	NEBULIZER COMPATIBILITY					Special pack available
			RF9	RF7 Dual Speed Plus	RF7 Dual Speed	RF8 Basic	RF8 Dual Speed	
 ACO611P	Mask soft Adult + elastic band (B022)	Single pack	✓	✓	✓	✓	✓	50 pieces package without elastic band (AC0473P)
 ACO612P	Mask soft Pediatric + elastic band (B023)	Single pack	✓	✓	✓	✓	✓	50 pieces package without elastic band (AC0472P)
 ACO437P	Mask soft Adult + elastic band (B022)	Single pack	✓	✓	✓	✓	✓	-
 ACO438P	Mask soft Pediatric + elastic band (B023)	Single pack	✓	✓	✓	✓	✓	-
 ACO462P	Mask soft Adult and Pediatric + elastic band (B022 + B023)	Single pack	✓	✓	✓	✓	✓	-

## Silicone Masks by Flaem - Baby




Full range of soft & safe silicone masks for babies. The masks can be connected with the universal L-bend connector, that is allowing to reach the baby in every position.

Code	Description	Package type	NEBULIZER COMPATIBILITY					Special pack available
			RF9	RF7 Dual Speed Plus	RF7 Dual Speed	RF8 Basic	RF8 Dual Speed	
 ACO389P	FLAEM BABY Mask Kit (Mask size 0 with bend) Size 0 = Neonatal - 1 month (B002 + B014)	Single pack with manual	✓	✓	✓	✓	✓	-
 ACO390P	FLAEM BABY Mask Kit (Mask size 1 with bend) Size 1 = 1-11 months (B002 + B015)	Single pack with manual	✓	✓	✓	✓	✓	-
 ACO391P	FLAEM BABY Mask Kit (Mask size 2 with bend) Size 2 = 1-3 years (B002 + B016)	Single pack with manual	✓	✓	✓	✓	✓	-
 ACO013B	L-Bend BABY connector (B002)	50 pieces	✓	✓	✓	✓	✓	-

## Mouthpieces

Flaem original mouthpieces	Code	Description	Package type	NEBULIZER COMPATIBILITY					Special pack available
				RF9	RF7 Dual Speed Plus	RF7 Dual Speed	RF8 Basic	RF8 Dual Speed	
 ACO441P	ACO441P	Mouthpiece with expiratory valves (B009)	50	✓	✓	-	-	-	-
 ACO427P	ACO427P	Mouthpiece (B017)	50	✓	-	✓	✓	✓	-
 ACO250P	ACO250P	Mouthpiece+Nosepiece Non-Invasive with valves (B009 + B010)	1	✓	✓	-	-	-	-

## Nosepieces

Flaem original mouthpieces	Code	Description	Package type	NEBULIZER COMPATIBILITY					Special pack available
				RF9	RF7 Dual Speed Plus	RF7 Dual Speed	RF8 Basic	RF8 Dual Speed	
 ACO448P	ACO448P	Nosepiece Non-Invasive ADULT/PEDIATRIC (B010)	50	✓	✓	-	-	-	-
 ACO450P	ACO450P	Nosepiece Standard ADULT (B026)	50	✓	✓	✓	✓	✓	-
 ACO451P	ACO451P	Nosepiece Standard PEDIATRIC (B027)	50	✓	✓	✓	✓	✓	-

## Filters - Tubing - Controls

Flaem original mouthpieces	Code	Description	Package type	NEBULIZER COMPATIBILITY					Special pack available
				RF9	RF7 Dual Speed Plus	RF7 Dual Speed	RF8 Basic	RF8 Dual Speed	
 ACO428P	ACO428P	Drug Saver Manual control	50	✓	✓	✓	✓	✓	-
 ACO225P	ACO225P	Nebulizer tube 2 mt. with fittings (B007)	1	✓	✓	✓	✓	✓	25 pieces (AC0457P)
 ACO35	ACO35	Nebulizer tube 1 mt. with fittings (B001)	1	✓	✓	✓	✓	✓	50 pieces (AC0432P)
 ACO578P	ACO578P	Nebulizer tube 2mt. with Long-therapy fittings. With TPU fitting (B033)	1	✓	✓	✓	✓	✓	25 pieces (AC0735P)
 ACO164P	ACO164P	Kit 2 Filters	2	-	-	-	-	-	50 pieces (AC0455P)
 ACO439P	ACO439P	Kit 2 Filters Basic	2	-	-	-	-	-	50 pieces (AC0485P)
 ACO45	ACO45	Kit 2 Filters	2	-	-	-	-	-	50 pieces (AC0456P)