

# Flaem Nebulizer systems



Portable / Silent



Intensive use



Frequent use



Seasonal use

NAME	Koala	LightNeb	Wi.Neb go	ProfiNeb
CODE	KO50P00	MF32E00	WI02P00	MI06P00
BUILT-IN STORAGE	No	No	No	No
AREA				
NEBULIZER				
COMPRESSOR	12V DC	Vibrating Mesh Technology	12V DC	12V DC
AIR OUTPUT	9 L/min'	-	9 L/min'	9 L/min'
USE	30' ON 30' OFF	5' ON - Automatic Switch Off	30' ON 30' OFF	30' ON 30' OFF
DIMENSIONS	W 16 cm H 17,7 cm D 13,5 cm	W 5 cm H 13,7 cm D 6,2 cm	W 13 cm H 5,5 cm D 9,5 cm	W 12 cm H 6 cm D 12 cm
WEIGHT	0,8 Kg	0,23 Kg (device) 0,12 Kg (batteries)	0,55 Kg (device) 0,14 Kg (batteries)	0,58 Kg

NAME	AirPro 3000 Plus	NebulFlaem 4.0	Hospineb	NebulAir+
CODE	AP50P00	Y301P00	FJ38P00	EL37P00
BUILT-IN STORAGE	No	Yes	No	Yes
AREA				
NEBULIZER				
COMPRESSOR	F3000	F2000	F2000	F2000
AIR OUTPUT	14 L/min'	14 L/min'	14 L/min'	14 L/min'
USE	Continuous	Continuous	Continuous	Continuous
DIMENSIONS	W 26 cm H 23,5 cm D 12 cm	W 20 cm H 10 cm D 30 cm	W 20 cm H 11 cm D 18 cm	W 20 cm H 10 cm D 30 cm
WEIGHT	2,4 Kg	2,4 Kg	1,8 Kg	2,4 Kg

NAME	RespirAir	Delphinus
CODE	FU02P00	DF30P00
BUILT-IN STORAGE	No	Yes
AREA		
NEBULIZER		
COMPRESSOR	F1000	F1000
AIR OUTPUT	10 L/min'	10 L/min'
USE	Continuous	Continuous
DIMENSIONS	W 16 cm H 10 cm D 16 cm	W 18 cm H 10 cm D 30 cm
WEIGHT	1,45 Kg	2,1 Kg

NAME	4Neb	AlphaNeb 4.0	Neb-Aid / Neb-Aid Duo	MasterNeb	Puppy	AirMate / AirMate Duo	AeroFlaem
CODE	CF10P00	AN11P00	JN33P00 / JN37P00*	PT13P01	KO00P00	SC36P00 / SC44P00*	NA15P00
BUILT-IN STORAGE	No	No	No	No	No	No	Yes
AREA							
NEBULIZER							
COMPRESSOR	F400	F400	F400	F400	F400	F400	F400
AIR OUTPUT	9 L/min'	9 L/min'	9 L/min'	9 L/min'	9 L/min'	9 L/min'	9 L/min'
USE	Continuous	Continuous	Continuous	Continuous	Continuous	Continuous	Continuous
DIMENSIONS	W 15 cm H 11 cm D 21 cm	W 18 cm H 13 cm D 19 cm	W 18 cm H 13 cm D 14 cm	W 16 cm H 11,5 cm D 16,5cm	W 16 cm H 17,7 cm D 13,5 cm	W 16 cm H 12 cm D 20cm	W 27 cm H 11 cm D 19 cm
WEIGHT	1,3 Kg	1,3 Kg	1,3 Kg	1,3 Kg	1,2 Kg	1,5 Kg	1,7 Kg

\* Duo version only