



# 03. SUC TION

## 03 Suction

Flaem Aspira is a range of high-quality professional vacuum systems that are ideally designed to deliver reliable suction, easy handling and cleaning. All models feature “Made in Italy” quality pumps, molding and high standard quality assembly. The housing with no gaps or grooves has been designed for easy cleaning and minimize risk of infections. The materials can withstand potent disinfection agents.



### Aspira Go

Choose the Aspira Go suction pumps by Flaem if you need a lightweight, reliable mobile airway suction pump for home healthcare. Available in single or double-pump version, Aspira Go models are powered by a removable Li-Ion battery pack.



### Aspira

Choose the Aspira suction pumps by Flaem if need lightweight, but robust and portable suction units. Available in three versions (11 lt, 13.5 lt, 20 lt) for different application needs.



### Aspira 40

Choose the new Aspira 40 for a new powerful and flexible solution. A new Made in Italy double-piston pump allows a vacuum flow up to 40 l/min for increased performance, fast suctioning, and high reliability. Aspira 40 has been designed to be used either transportable or rack version, because of the exclusive design of a trolley connection plate.



#### Heavy Duty

High quality 100% Made in Italy pumps, molding and high standard quality assembly.



#### Seamless design

Housing designed for easy cleaning and minimise risk of infections



#### Patented Jar-Holders

Patented high quality jar holders for a better portability and no vibrations.



#### Full control of vacuum level

Patented system for a precise control of the suction level.

# Aspira Go 30

Professional battery-operated portable suction device

Code AR30P00 | Model P1611EM-30 | Pack 2 pcs



Compact, lightweight battery-operated portable surgical aspirator. Easily transportable thanks to lithium batteries easily replaceable. Simple maintenance and easy sanitization. Easy and intuitive reading of the functions. Practical manual control of the suction flow.

# Aspira Go 20

Professional battery-operated portable suction device

Code AR20P00 | Model P1611EM-20 | Pack 2 pcs



Compact, lightweight battery-operated portable surgical aspirator. Easily transportable thanks to lithium batteries easily replaceable. Simple maintenance and easy sanitization. Easy and intuitive reading of the functions. Practical manual control of the suction flow.

Portable

- Made in Italy
- Adjustable vacuum level
- Runtime 30' ON - 30' OFF
- Silent
- Sterilizable jar and tubes
- Real flow 30 L/min'
- Jar-holder (1 included)



MD



Portable

- Made in Italy
- Adjustable vacuum level
- Runtime 30' ON - 30' OFF
- Silent
- Sterilizable jar and tubes
- Real flow 20 L/min'
- Jar-holder (1 included)



MD



## Accessories included



## Certifications



## Accessories included



## Certifications



# Aspira 20

Professional desk type suction device

Code AS20P00 | Model P1211EM/20 I | Pack 2 pcs



Designed for medical and surgical use in hospitals, clinics, and surgeries, this device is also ideal for home pharyngeal aspiration therapies. Its durable polymer construction ensures both robustness and easy portability. The snap-on jar holder accommodates either a single-patient 800 ml or a sterilizable 1000 ml graduated jar. Silicone tubes connect the jar to the suction device via a high-efficiency filter to prevent contamination.

# Aspira 13.5

Professional desk type suction device

Code AS01P00 | Model P1211EM/13,5 I | Pack 2 pcs



Designed for medical and surgical use in hospitals, clinics, and surgeries, this device is also ideal for home pharyngeal aspiration therapies. Its durable polymer construction ensures both robustness and easy portability. The snap-on jar holder accommodates either a single-patient 800 ml or a sterilizable 1000 ml graduated jar. Silicone tubes connect the jar to the suction device via a high-efficiency filter to prevent contamination.

## Desk Type

- Made in Italy
- Adjustable vacuum level
- Heavy duty
- Silent
- Sterilizable jar and tubes
- Real flow 20 L/min'
- Jar-holder (1 included)



## Accessories included



## Certifications



## Desk Type

- Made in Italy
- Adjustable vacuum level
- Heavy duty
- Silent
- Sterilizable jar and tubes
- Real flow 13.5 L/min'
- Jar-holder (1 included)



## Accessories included



## Certifications



# Aspira 11 Professional desk type suction device

Code AS01P00 | Model P1211EM / 11 | Pack 2 pcs



Designed for medical and surgical use in hospitals, clinics, and surgeries, this device is also ideal for home pharyngeal aspiration therapies. Its durable polymer construction ensures both robustness and easy portability. The snap-on jar holder accommodates either a single-patient 800 ml or a sterilizable 1000 ml graduated jar. Silicone tubes connect the jar to the suction device via a high-efficiency filter to prevent contamination.

## Desk Type

- Made in Italy
- Adjustable vacuum level
- Heavy duty
- Silent
- Single patient jar and tubes
- Real flow 11 L/min'
- Jar-holder (1 included)



## Accessories included



## Certifications



# Aspira 40

Flexible high flow suction device,  
portable or rack version

Code AS40P00 | Model P0720EM / 40 I | Pack Single Unit

COMING  
2025



Aspira 40 is a professional surgical aspirator intended for use in healthcare facilities such as hospitals, clinics and surgeries. It is also indicated for home therapies for pharyngeal aspirations. It can be installed on an optional trolley and equipped with a comfortable pedal for intermittent use. Flaem jars are equipped with float protection and a high efficiency hydrophobic filter for anti-microbial and antibacterial contamination. It comes with a suction canister on the left side, but it's easy to add an optional jar holder on the opposite side.

Flexible

Validato dal TÜV  
Flow  
40 l/min  
Validated by TÜV

- Made in Italy
- Heavy duty
- Double compressor allows to reach a flow of 40 L/min'
- Foot switch (optional) for intermittent operation
- Adjustable vacuum with gauge display
- Mode switch for continuous or intermittent use
- Useful trolley (optional) with connection plate for hosting liner containers

MD



## Accessories included







## Certifications

CE PENDING

PHYTHALATES & BPA FREE



## Optional Accessories for Aspira 40

	CODE	DESCRIPTION
	ACO772P	Foot switch for intermittent operation
	ACO773P	Trolley kit with connection plate (4 support rings included)
	15224A0	Additional jar holder
	ACO001A	Single-patient jar 800 ml



Aspira 40 can be installed on a trolley (optional) equipped with 5 spokes and swivel wheels, with safety brake. The attachment takes place via a special plate which is inserted into the bodies and is fixed in a stable position by two safety screws. The trolley (optional) is equipped with a special hub to fix up to 4 holder rings, useful for hosting liner containers (high vacuum/high flow).