



04.

DIAGNOSTICS & WELLBEING

04

Diagnostics & Wellbeing

The Diagnostics & Wellbeing section of the Flaem catalog showcases a range of products designed to empower individuals to take control of their health and well-being. These products support accurate diagnosis and promote a healthier lifestyle.



Thermometers

Flaem provides both manual and digital thermometers for precise temperature measurement. These thermometers are designed with user-friendliness and accuracy in mind, aiding in the early detection of potential health issues



Blood Pressure Monitors

Flaem's digital blood pressure monitors offer convenient and reliable blood pressure tracking. These monitors contribute to the management of cardiovascular health by providing essential data for individuals and healthcare professionals



Humidifiers

Flaem's hot-steam humidifiers help maintain optimal humidity levels in indoor environments. This can be particularly beneficial for respiratory health, especially during dry seasons or in arid climates.

These products collectively support individuals in monitoring their health, making informed decisions, and or taking proactive steps to improve their overall well-being.

Accu-Therm® View

Code TRM10 | Model CR.W00 | Pack 10 pcs



Clinical manual thermometer, completely disinfectable thanks to its glass materials. A shock-free case protects the thermometer when it is unused, while easing the temperature reset by a swinging movement. Environment friendly: no batteries needed and non-toxic mercury-free liquid (galinstan).

Accu-Therm®

Hard Tip → Code TRM07 | Model DMT-2031 | Pack 12 pcs
Flex Tip → Code TRM12 | Model DMT-4726 | Pack 12 pcs



Flaem Accu-Therm® digital thermometers use micro-computerized technologies that guarantee the quality of the products. The maximum accuracy of the measurements is guaranteed by a test to verify that these instruments perform at each switch.

Accu-Therm® Hard Tip

- Temperature detection in 60"
- LCD display
- Memory from the last measurement
- Acoustic signal end-measurement
- Automatic power-off
- Low battery signal

Accu-Therm® Flex Tip

- Digital thermometer with flexible tip
- Temperature detection in less than 30"
- Hypoallergenic golden soft tip "Nikel free"
- Memory from the last measurement
- Wide display with colour switch (red: fever – green: ok)
- Acoustic signal end-measurement at a temperature of 37.8°C
- Low battery signal

- Ca Thermometer with Galinstan (mercury-free)
- Ni/PVC Hypoallergenic, Nickel and PVC free
- ✓ Resistant shock-free case
- Extremely easiness of use (no-battery operated)
- Italian design



- 30 Detection of temperature in less than 30 seconds model TRM12
- 60 Detection of temperature in 60 seconds model TRM07
- Acoustic signal at end measurement
- Ni Hypoallergenic, Nickel Free
- Auto power off



Packaging



Manufacturer

Made in PRC
Wuxi Medical Instrument Factory
Address: No. 86 Dongjie, Zhangjing, Xibei Town,
Wuxi City, Jiangsu, 214194 China.

CE 1023 EC REP Lotus Global Co., Ltd
Address: 15 Alexandra Road, London
UK nw80DP

Packaging



Manufacturer

Made in PRC
JOYTECH HEALTHCARE Co., Ltd.
No. 365, Wuzhou Road, Yuhang
Economic Development Zone,
Hangzhou city, 311100 Zhejiang, China

CE 0197 EC REP Shanghai International
Holding Corp. GmbH
(Europe) Eiffestrasse 80,
20537 Hamburg, Germany

SfigmoLife Pro

Code SFG62 | Model TMB-1490-CS | Pack 12 pcs



Digital blood pressure monitor with MWI technology. Compact and user-friendly, it offers rapid readings thanks to MWI technology. Equipped with a large backlit display for optimal readability, it can store up to 250 measurements per user. Includes arrhythmia detection and cardiac risk classification according to WHO guidelines. Ideal for regular home monitoring.







SfigmoLife Evo

Code SFG61 | Model TMB-1491-S | Pack 12 pcs









Digital blood pressure monitor, leverages MWI technology for rapid, accurate readings. Its ESH-certified quality ensures reliability. With dual-user memory (60 readings each), arrhythmia detection, average measurement tracking, and a watch mode, the SfigmoLife EVO is a versatile home health monitoring tool.

 ESH validation

-  2 memory channels up to 500 measurements
-  4 AAA batteries included
-  Cardiac-risk classifications
-  Irregular heartbeat indicator
-  Travel bag included
-  Backlight display



-  Clinically tested for the efficacy
-  2 memory channels up to 120 measurements
-  4 AAA batteries included
-  cardiac-risk classifications
-  Irregular heartbeat indicator
-  Travel bag included



Manufacturer

Made in PRC
 0123
 EC REP
 Guangdong Transtek Medical Electronics Co., Ltd.
 Zone A, No.105, Dongli Road,
 Torch Development District, Zhongshan,
 528437, Guangdong, China.



M-Sized cuff (22 cm – 32 cm): within accessories
 L-sized cuff (22 cm – 42 cm): purchasable - code 18312

Manufacturer

Made in PRC
 0123
 EC REP
 Guangdong Transtek Medical Electronics Co., Ltd.
 Zone A, No.105, Dongli Road,
 Torch Development District, Zhongshan,
 528437, Guangdong, China.



M-Sized cuff (22 cm – 32 cm): within accessories
 L-sized cuff (22 cm – 42 cm): purchasable - code 18312

