

Smarty
 Compact, quite,
 hand-held nebulizer
 Made in Italy



**Aerosol on demand
 or continuous mode**



**Comfortable
 adult and paediatric
 Soft Touch masks**



**Power supply: 2 AA
 batteries, Alkaline or
 Rechargeable**



**The respiratory
 therapy**

— 50 —
 years of italian excellence

Smarty



Compact, quite,
hand-held nebulizer
Made in Italy



Technical Features

VMT Vibrating Mesh Technology nebulizer.

It nebulizes also suspension drugs such as beclomethasone dipropionate and budesonide Continuous.

Power supply/Power	2.4 - 3VDC - 1.5W 2 AA batteries, alkaline or rechargeable	
Maximum fill volume (l)	8 ml	
Transducer frequency	172 KHz	
Sound level (at 1 m)	Less than 20 dBA	
Device dimensions	5.4 (length) x 6.9 (width) x 14 (height) cm 2.13 (length) x 2.72 (width) x 5.5 (height) in	
Device Weight	130 g - 4.58 oz	
Carrying bag dimensions	5.9(length) x 3.15(width) x 5.9(height) in	
Total weight (device, carrying bag, batteries and accessories)	330 g - 11.64 oz	
LED	Blue: ON - Blinking red: low batteries	
Automatic switch off	After 5 minutes or in case of accidental opening or no medication present	
Approvals	In compliance with 93/42/EEC Medical Device Directive	 
	Complies with: European standard EN 10993-1 "Biological Evaluation of Medical Devices" and European Directive 93/42/EEC "Medical Devices" Phthalates free. In conformity with Reg. (EC) no. 1907/2006	

Nebulisation technical features

Delivery ml/min	≥ 0.25 (1)
Characterization *MMAD	4.00 µm approx (2)
Breathable fraction < 5 µm (FPF)	65% approx (2)

(1) The nebulising speed has been measured with a 0.9% saline solution at 23° C according to the internal Flaem procedure I29-P07.5. It may vary according to the nebulising head supplied, medicinal product and environmental conditions. Values shown for the aerosol delivery may also vary depending on the patients respiratory Capacity.

(2) In vitro testing performed by University of Parma in compliance with the EN 13544-1:2007 + A1 European Standard for aerosol therapy devices Standard. More details are available on request.

Accessories

Soft silicone mouthpiece
Soft Touch Adult mask
Soft Touch Child mask
2 AA alkaline batteries
Convenient carrying bag

Smarty

Code	MF02E00
Model	P031SEM
Master carton	8 pieces



To know more



1701600 © Copyright by FLAEM NUOVA S.p.A., 2020 - FLAEM NUOVA reserves the right to make technical and functional alterations to this product without prior notice.