



# 01. NEBU LIZA TION

- 01.A Rhino Clear Range
- 01.B Nebulizer Systems
- 01.C Nebulizers & Kits
- 01.D Accessories

## 01.A Rhino Clear® Range

A complete range of innovative devices, standalone or to be connected to a compressor nebulizer, specifically designed for a correct treatment of the upper respiratory tract.

The Rhino Clear® range is:



### Patented, 2 chamber system

Different chambers for drug and saline solution, and for secretions collection.



### Comfortable and safe

Ergonomic design, soft adapters for adult and children noses.



### Effective

Nebulizes 10 ml in just 2 minutes and removes mucus, dust and harmful substances.



### Portable

Rhino Clear Sprint and Rhino Clear Mobile are designed for use on the go.



# Rhino Clear® / Rhino Clear® Ergo Plus

- Rhino Clear® | Model RC1-1**
- Neutral Code ACO556P
  - Blu Code ACO415P
  - Yellow Code ACO416P

**Rhino Clear Ergo Plus | Model RC2-2**  
Code ACO702P



Rhino Clear® is a fast and effective system to clean the upper respiratory tract. Paired with a compressor nebulizer, it removes secretions, dust and harmful substances that have accumulated inside the paranasal sinuses and collects them in a special chamber. It can be used with saline solutions (isotonic, hypertonic, liposome solutions, etc.), thermal waters or medications commonly available in pharmacies.

# Rhino Clear Sprint

Code DR00P00 | Model P0709EM | Pack 6 pcs



Portable nasal wash, with batteries. Rhino Clear Sprint is Effective (correct particle size for nebulization in the upper airways), Comfortable (non-invasive, easy to clean, can be used everywhere) and suitable for both adults and children. Can be used with saline solutions (isotonic, hypertonic, liposome solutions, etc.), thermal waters or medications commonly available in pharmacies.

**Made in Italy** MD

- Compatible with the entire range of Flaem compressor nebulizer systems
- Start/Stop button
- Proven effectiveness
- Fast treatment (2 minutes for 10ml nebulization)



**Made in Italy** MD

- 2 AA alkaline batteries included (or rechargeable)
- Nasal soft-touch adaptors for adults and children
- Low energy consumption
- Fast treatment (2 minutes for 10ml nebulization)



## Rhino Clear Ergo Plus Accessories included



Up to 3 years      3-12 years      12+ years

## Rhino Clear® and Rhino Clear Ergo Plus technical features

MMD µm (Mass Median Diameter) <sup>(1)</sup>	> 10
Delivery ml/min <sup>(2)</sup>	5
Maximum fill volume	10 ml

(1) Values detected with the Malvern Spraytec instrument and referring to the use of the nasal wash with 0.9% NaCl physiological solution. Particle sizes greater than 10 µm are suitable for the treatment of upper airways and in particular for cleaning the nasal cavities. - (2) Average values found with 0.9% NaCl saline solution according to internal Flaem procedure.

## Rhino Clear Accessories included



Up to 3 years      3-12 years      12+ years

## Certifications



## Rhino Clear Sprint technical features

MMD µm (Mass Median Diameter) <sup>(1)</sup>	> 10
Delivery ml/min <sup>(2)</sup>	5
Maximum fill volume	15 ml

(1) Values detected with the Malvern Spraytec instrument and referring to the use of the nasal wash with 0.9% NaCl physiological solution. Particle sizes greater than 10 µm are suitable for the treatment of upper airways and in particular for cleaning the nasal cavities. - (2) Average values found with 0.9% NaCl saline solution according to internal Flaem procedure.

## Accessories included



Up to 3 years      3-12 years      12+ years

## Certifications



# Rhino Clear Mobile

Code DR04P00 | Model P0907EM | Pack 8 pcs



Induction charging station

Portable nasal wash, with internal micro-compressor and rechargeable battery. Rhino Clear Mobile is Effective (correct particle size for nebulization in the upper airways), Comfortable (non-invasive, easy to clean, can be used everywhere) and suitable for both adults and children. Can be used with saline solutions (isotonic, hypertonic, liposome solutions, etc.), thermal waters or medications commonly available in pharmacies.

Made in Italy

MD

- Battery Ni-Mh low discharge with charging station
- Nasal soft-touch adaptors for Adult and Child
- Low energy consumption
- Fast treatment (2 minutes for 10ml nebulization)
- Travel bag included



## Rhino Clear Mobile technical features

MMD $\mu\text{m}$ (Mass Median Diameter) <sup>(1)</sup>	> 10
Delivery ml/min <sup>(2)</sup>	5
Maximum fill volume	15 ml

(1) Values detected with the Malvern Spraytec instrument and referring to the use of the nasal wash with 0.9% NaCl physiological solution. Particle sizes greater than 10  $\mu\text{m}$  are suitable for the treatment of upper airways and in particular for cleaning the nasal cavities. - (2) Average values found with 0.9% NaCl saline solution according to internal Flaem procedure.

## Accessories included



Up to 3 years    3-12 years    12+ years

## Certifications

CE 0051



PHYTHALATES & BPA FREE

# Rhino Clear Aspirator

Code ACO423P | Model J002 | Pack 1 pcs · Code ACO510P | Pack 1x6 pcs



The Rhino Clear aspirator gently removes mucus from your child's nose with a simple push and release. Its user-friendly design and safe nozzle make it easy for parents to use. It is also easy to clean and sanitize.

Made in Italy

MD

- Patented mouthless suction
- Simple operation
- Gentle on the nose
- Easy to clean



## Packaging



## Accessories included



## Certifications

CE

PHYTHALATES & BPA FREE